


Fix-GPC UV Chromatogram of the Calibration Solution

Printing time: Wed Dec 18 17:06:08 2013
Injection time: Wed Dec 18 15:48:48 2013

File opened: D:\Documents and Settings\Yan\Desktop\GPC standard test
datas\Hongkong GPC EPA 3640 Test Report - R1.hw
Multiplying factor: 1, Dividing factor: 1

Baseline noise: 11637uV, baseline drift: 584947uV



Rank	Time	Name	Area	Resolut.	Corre. f
1	6.475	corn oil	35475206	1.60	
2	9.340	Bis(2-ethylhexyl)phthalate	22685790	1.06	
3	11.302	Methoxychlor	21080318	2.03	
4	15.605	Perylene	19729114	1.01	
5	17.995	Sulfur	13710566	0.00	
Total			112680994		